



The ERWEKA Disso.NET 4 is the perfect 21 CFR Part 11 compliant companion to all our dissolution systems, ranging from dissolution offline systems over online-to-on-/offline systems, up to the fully automated RoboDis II+ system.

The software helps you with standard dissolution jobs, easily handles qualification tasks and provides control over each single function and connected devices (e.g. dissolution tester with UV-Vis spectrophotometer). After finishing those tests, Disso.NET creates extensive reports with corporate logos for PDF-file export or exports your results in XLS or XML.


Our extensive Audit Trail according to latest 21 CFR part 11 thoroughly documents all changes done to the system (what, who, when and why) and can be easily searched and filtered by the Audit Trail Viewer.

Disso.NET 4 has an active directory connection. This allows cross-system login data and user passwords, thus avoiding countless passwords and ensuring easy handling.



FROM OUR DISSOLUTION TESTERS TO THE ROBODIS II FULL CONTROL FOR ALL OUR SYSTEMS

Disso.NET 4 allows control over all our dissolution systems. Simpler dissolution online systems, complex dissolution on-/offline with HPLC Systems and the complex RoboDis II+ can all be controlled by Disso.NET 4, while reducing device and system complexity to a few intuitive screens with a logical top-to-bottom approach.



[click to enlarge](#)

DISSOLUTION TESTER SOFTWARE DISSO.NET 4 KEY FEATURES



EASY DATA MANAGEMENT

Users, methods, reports and more can be easily managed with Disso.NET 4. Multi-user support with different user roles (i.e. Operator, Senior Operator, Administrator - all configurable) allows for fine access control to the system and data and can be even fine-tuned with single permissions. The automatic log-off rounds up the access safety features and makes sure that only approved personnel has access.

SYSTEM SUITABILITY & MULTITASKING

Constant system suitability tests ensure highest reliability by constantly checking the system. The JobController principle allows easy multitasking while testing.

DATABASE WITH BACK-UP AND RESTORE

DOCUMENTATION & REPORTS

Disso.NET 4 fully documents tests: With the click of a button, customers can easily create reports with customized logos - from compact analysis reports to method data, calibration data or raw data. The integrated workflow enables easy approval and review before publishing reports for the paperless lab.

LANGUAGE SUPPORT

Language support for English and German is already built into Disso.NET 4 dissolution tester software. More languages can be easily added via language files.

POWERFUL MODULES

DISSO.NET EDITOR

All main features of Disso.NET 4 are located in separate windows, ensuring a clear focus on the task at hand. These editors can be used to change media, edit filters and create new measurement, evaluation and medium preparation methods.

Even for UV calculations, Disso.NET offers a unique interface specifically focused on its main task.

